SPENQ • MEDICAL

COMPANY PROFILE





Anji SPENQ Industrial Co., Ltd is a professional export-oriented enterprise that major in medical dressing and other medical products. Our group is formed through the amalgamation of 4 factories and a marketing centre. It is located in Anji, in China's Zhejiang province, which is famously known as "Hometown of Medical Dressing in China".

We are proud of our unrivaled product range and our ability to rapidly customize our products to cater for each customer's individual needs. Our group roundly actualizes the international quality managing system of ISO13485:2003, and got CE and FDA certificate.

Our values are: excellence, professionalism, and innovation. When translated into practice, it means that we strive to achieve excellent product quality, swift service both before and after sales, specialized product knowledge, constant innovation, and organizational amalgamation.

At SPENQ, we believe the main idea "customer central". Our clients are the foundation of our business, and their satisfaction is the only measure of our success. It is for this reason that we established an unique customer service department in order to expediently contact with our customers.

At SPENQ, our mottos are "communication leads to development" and "cooperation leads to excellence". We will like to extend our welcome to anyone from all walks of life, from both at home and abroad, to join in our adventure.





PRODUCTION

SPENQ Medical Owned:

- 7 Professional Assembly Lines
- 380 Sets Automatic Machines
- 32000M² Factory Area
- 35 Experienced Engineers1200 Skilled Workers



CERTIFICATION



















View more products at:

<u>WWW.SPENQ.COM</u>







MAD1 High Elastic Bandage o







Polyester + Nature latex



Skin White



5cm x 4.5m	
7.5cm x 4.5m	
10cm x 4.5m	
15cm x 4.5m	



Characteristics:

- O High intensity and stretch, the extensibility is over 295%.
- O Easy using.
- O Washable. Description 🖸

Sports, emergency saving, training, military ambulance, air tragedy rescue and so on. Polyester/rubber knitted construction, ideal for high compression treatment of condition.



MAD2 Elastic Crepe Bandageo







Cotton + Spandex



Nature Bleached Skin



5cm x 4.5m 7.5cm x 4.5m 10cm x 4.5m 15cm x 4.5m



Characteristics:

- O High elasticity, extensibility is over 180%.
- O Easy using.
- O Washable.

Description 🕽

Applications:
For treatment, after-care and prevention of recurrence of working and sports injuries, after-care of varicose veins damage and operation as well as for therapy of vein insufficiency.

SPEND

BANDAGE



MAD3 Cotton Crepe Bandage o



Material 0	Type 1: 100% cotton. Type 2: 98% cotton + 2% spandex.
Color	CNature C Bleached
(Tames)	5cm x 4m
FI 0	7.5cm x 4m
Size D	10cm x 4m
	15cm x 4m



Characteristics:

- Wonderful feeling, comfortable to wear.
- Breathing, very good absorbency
- Washable.
- Description 🔾

Applications:
For all of after-care of operation, injuries, soft tissue fractures and joint swelling, etc.

MAD4 Elastic Plain Bandage o







OBleached OSkin Nature



5cm x 4.5m 7.5cm x 4.5m 10cm x 4.5m 15cm x 4.5m



Characteristics:

- O High elasticity, extensibility is over 200%.
- O Easy using.
- Washable.

Description 🔾

Applications:

The elastic plain bandage is specially designed for protecting skin wounds and changing wound dressings. It is extremely comfortable to wear and convenient to use.

Conforming Bandage (Thin PBT)





100% Polyester



Bleached



2.5cm x 4m 10cm x 4m 15cm x 4m 5cm x 4m 7.5cm x 4m



Characteristics:

 Soft and comfortable, with high elasticity. The joint will not be confined after use PBT, will not cause bad blood circulation and movement of the joint.

Description 0

Applications: Ideal for binding all parts of the body, outdoors


practice, first-aid swaddle and soft tissue damage.





MAD5 Ideal Bandage (Thick PBT)





9


Description 0

Soft and comfortable, with good elasticity, better intensity than confirming bandage.
 The joint will not be confined after use PBT, will not cause bad blood circulation and movement of the joint.

Applications:

Ideal for binding all parts of the body, outdoors practice, first-aid swaddle and soft tissue damage.

MAD7 Gauze Bandage with cutting edge (W.O.W Bandage)



100% cotton (40's/40's) Material 0 Color D Bleached 13threads/19x15mesh 17threads/26x18mesh Mesh 0 20threads/30x20mesh 2.5cm x 5m 10cm x 5m 5cm x 5m 15cm x 5m Size 0 7.5cm x 5m Characteristics: Soft and comfortable. O Absorbent and breathing. Description 0





Material 0	100% cotton (40's/40)'s)
Color O	○ Bleached	
<u>re</u>	13threads	/19x15mesh
Table 1	17threads/26x18mesh	
Mesh	20threads	/30x20mesh
	2.5cm x 5m	10cm x 5m
	5cm x 5m	15cm x 5m
Size D	7.5cm x 5m	
Description 0	Characteristics: Soft and comfortable Absorbent and breath Woven edge knitted co	

better intensity.

SPEND

BANDAGE



MAD9 Plaster of Paris Bandage o







Gauze + Plaster



White



Description 0

8	5 cm		2.7m
	7.5cm	5	78130
Width	10 cm	E	3 m
_	15cm	5	3.6m
	20cm		3,0111



Characteristics:

- Solidification time is fast and steady: 2-5 minutes.
- Time of soaking in water is fast: 5-10 seconds.
- It becomes dry very fast after solidification. The intensity reaches 3.2Mpa after 8min.
- Its caloric value is compatible with human temperature: 35-40c.

Applications:

It is used for fracture orthopedic and modelmaking. It is safe and easy to be used and has won the favour from doctors and patients.









Type 1: Fiberglass Fabric + High Polymer Glue Type 2: Polyester Fabric + High Polymer Glue



○ White Blue Green Pink CYellow Orange Pattern



2.5cm x 1.8m 10 cm x 3.6m 5 cm x 3.6m 12.5cm x 3.6m 7.5cm x 3.6m 15 cm x 3.6m



Description 🔾

Characteristics:


- O High intensity(20 times of P.O.P bandage)
 O Very light(one fifth of P.O.P bandage)
- Good breathing, can avoid canker of the skin, and washable
- The penetration rate for X-ray is 100%.

It is a new orthopedic bandage instead of P.O.P bandage. With advanced technology, it enjoys good quality at the reasonable price, which are widely accepted by orthopedists in

4













Description 0

Characteristics:

 High intensity, and has good shaping and stability.

O Easy to use.

O Convenient to observe and adjust.

It is one of new surgical orthopedic products. Ideal for surgery, fixing parts, and model-



MA12 Orthopedic Padding o





Material 0

Type 1: 100% Polyester.

Type 2: 100% Cotton.

Type 3: 100% Viscose.

Type 4: 50%viscose + 50%polyester.

Color 0

White



5 c m 7.5cm 10 cm 15cm 20cm

2.7m Length 3 m 3.6m



Characteristics:

Soft, smooth and conformable to the body's shape.

Description 0

Applications: Orthopedic Padding is the wide-used supporting material for surgery and swaddling of P.O.P bandages and orthopedic casting tapes. It can be used under all types of casting material, providing effective, comfortable and safe casting.



Tubular Bandage (with latex)







Cotton + Rubber



Size 0

Nature



- #A: 4.5cm infant feet or arms.
- #B: 6.5cm small hands or limbs.
- #C: 6.75cm adult hands, arms or legs.
- #D: 7.5cm large arms or legs.
- #E: 8.75cm legs or small thighs.
- #F: 10cm large knees or thighs.
- #G: 12cm large thighs.
- #H: 15cm extra large thighs.
- #I: 17.5cm small trunks or extra large thighs.
- #J: 21.5cm medium trunks.
- #K: 26cm medium large trunks.
- #L: 32.5cm large trunks.
- Note: We can do other sizes according to customers' req.



Characteristics:

- Breathable and comfortable.
- O Highly absorbent and well tolerated.
- Washable.

Applications: It is designed as a pre-wrap under a cast, and give ideal pressure and support over the area.



Stockinet Bandage (without latex)







100% cotton.



Color 0

Size 0

Nature



- #A: 5cm infant feet, arms or small hands.
 - #B: 7.5cm adult hands, arms or legs.
- #C: 10cm legs, large knees or thighs.
- #D: 15cm large thighs.
- #E: 20cm extra large thighs or trunks.
- Note: We can do other sizes according to customers' requirements



Characteristics:

- Breathable and comfortable, no latex.
- Highly absorbent and well tolerated.
- Washable.

Description 🗘

Applications:

It is designed as a pre-wrap under a cast, and give ideal pressure and support over the area.

Description 0

SPENI







MA15 Tubular Net Bandage o







Polyester + Rayon + Nature latex



White



0.8cm - finger or toe.

1.7cm - hand, wrist, ankle or lower arm.

2.0cm - hand, elbow or upper arm.

#3: 2.2cm - hand or foot.

2.6cm - foot or arm.

3.0cm - foot, arm or knee. 3.8cm - arm or leg.

#5:

Size 0

#5.8: 4.4cm - leg or head of children.

5.5cm - head of adult, shoulder, oxter, hip or venter,

#7: 6.8cm - shoulder, oxter, back, hip or venter.

8.0cm - shoulder, oxter, back, hip or venter.

For elastic hat, we provide sizes: #5.0, #5.5, #5.8

Note: We can do other sizes according to customers' requirem



Characteristics:

Economic: Few material but cover large area.
 Breathing, for a fast recovery.

Match any part of the body, fit well.

for the place that general bandage can't match.



MA16 Cohesive Elastic Bandage o







Type 1: Non-woven fabric + Nature latex Type 2: Cotton fabric + Nature latex

Color C

○ White ■ Blue ■ Green @ Pink ○ Yellow ●Orange ●Black ●Red ○Pattern



Description 🔾

2.5cm x 4.5m 10cm x 4.5m 5cm x 4.5m 15cm x 4.5m

7.5cm x 4.5m



Characteristics:

O High elasticity, self-adhering, not stick to hair, skin, clothing, no clips needed.

Easy hand tear, no scissors needed.

Soft, breathable and comfortable.

O Waterproof. No slip when meet sweat or water.

Applications:

Widely used in medical and sports area, also for veterinary and pet wrap. It can be used for fixing basic dressings or supporting muscle avoid sprain and strain.




Description 🖸

Quickly, simple, and easy for anyone, anywhere, anytime.

Fit all parts of the body immobilization, especially

MA17 First Aid Bandage o







Type1: Conforming bandage + Absorbent pad Type2: Conforming bandage + Ties



Bleached



4cm x 6cm 6cm x 8cm 8cm x 10cm 10cm x 12cm



Characteristics:

 High absorbent and non-stick to skin. Easy to take care of the dressing from the wound.

Description 0

Applications:

It is used when first treatment of middle size skin, wound and cutting injuries, it can be supplied according to customer specification.



MA18 Triangular Bandageo







Type 1: 30/40 g/m2 non-woven. Type 2: 40's cotton.



Color 0

C Bleached

ONature



Size 0

90x90x127cm 96x96x135cm



Characteristics:

O Very light and compact. With safety pins is available.

Description 0

Applications: It's designed to hang the arm, keep splints in place, or enswathe the wound in head. Mainly used on injured hand or arm.

SPENO















9

Description 🔾

Characteristics:

Low breathing resistance with high filtration

 Easy-breathing exhalation valve with inner silicon-piece can facilitate air exchange and increase people comfort.

Applications: It is widely used for surgical operation, fastening dressing and catheter, sports protection etc.



MB02 Non-woven Tape o

(Microporous Tape)









Characteristics:

- Soft and comfortable, hypoallergenic.
- Good adhesiveness to skin.
- Easy to tear, making usage very convenient.

Description 🔾

Applications:
A gentle, breathable tape for sensitive skin, recommended especially for pediatric and geriatric patients. Used for secures dressings, tubes, lines and I.V. catheters.











PE film + Hot-melt glue



Transparent













Description 🔾







Characteristics:

- Soft and comfortable, hypoallergenic.
- Low sensitive, the adhesive ability is good.
- Transparent and clean.

Applications:
A gentle, breathable tape for sensitive skin, recommended especially for pediatric and geriatric patients. Used for secures dressings, tubes, lines and I.V. catheters.










Silk-like Fabric + Hot-melt glue




O White









Len	5	Υ
gth	10	Υ



Characteristics:

- Soft and comfortable, hypoallergenic.
- Sawtooth edge, making it easy to tear.
- Good permeable and Waterproof.

Description 0

Size 0

Applications:

Ideal for bulky dressings that require high conformability, also for post-operative surgical dressings.

SPEND














Cotton fabric with polyethylene coating + Hot-melt glue




White







	1.25	cm	
Width	2.5	cm	E
	5	cm	P
	7.5	cm	5
	10	cm	
	10	CIII	





Characteristics:

Waterproof, oil resistant and easy to write on.
 High strength, and will not permit moisture to penetrate into skin or dressing.

Hypoallergenic and latex-free.

Description 0

Applications: It is widely used for surgical fastening dressing and low-weight catheter or as the fastening of nursing wound after large-area operation.



MBD6 Sport Tape o







(3)	8	2.5 cm	E	10Y
Size 0	Width	3.75cm	png	10m
- SHELL	_	5 cm	5	15Y



Characteristics:

Air permeability and water-resistant.High strength.

Sawtooth edge, making it easy to tear.

Description 🗘

Applications:

It is designed for medical, health care, sports and other areas.



10 MB07 Elastic Adhesive Plaster o







Elastic cotton fabric + Hot-melt glue.



White C Yellow Blue Red Skin Black



2.5cm x 4.5m 5cm x 4.5m 7.5cm x 4.5m 10cm x 4.5m



Characteristics:

- Good elasticity.
- Good adhesiveness, hypoallergenic.
- O High tensile strength.

Description 🐧 Applications:

Ideal for dynamic/pressure dressings requiring elasticity, support strapping, fractures and muscle injuries. Also can control and reduce post-injury edema by providing firm support and compression.

10

MBDB Elastic Adhesive Bandage o







100% Cotton Cloth + Zinc Oxide Glue



Color 0

White Skin



Size 0

Description 🗘

5cm x 4.5m 7.5cm x 4.5m 10cm x 4.5m



Characteristics:

- O Porous adhesive to allow skin to breathe for maximum adherence.
- Central yellow line for correct alignment of bandage to anatomy.
- O Non-fray fluffy-edged cloth that prevents cutting for maximum comfort (white products only).

Applications:
Provide compression support for joints during high-stress sport activity. The bandages help minimise swelling and provide additional support as an overwrap. Suitable for compression bandaging of ankles, knees, fingers, shoulders and elbows.








Spun-laced Wound Dressing Tape o





Material O	Spun-laced non-woven + Hot-melt glue.

Color O	○ White

Size O		5cm x 10m	
		10cm x 10m	
	Size O	15cm x 10m	
		20cm x 10m	

2		teristics: ble and comfortable.	
		20cm x 10m	
Size	ize O	15cm x 10m	
		Tocili x Tolli	

Description 0

Applications:
Convenience on areas with edema and on moving areas, especially for all kinds of surgical attempt areas, on shoulders, elbows, armpits, process areas, in post-operative usage.

Good adhesiveness, hypoallergenic.








Material 0	Polyurethane Foam
Color O	○ White ○ Skin ● Blue ● Red
Size 0	7cm x 30Y
Description •	Characteristics: O Soft and breathable O Good tensile strength O High elasticity Applications: It is used under adhesive sports tapes to help prevent irritation from repeated taping. It also aids in painless tape removal

WOUND DRESSING

MED1 Wound Plaster O







Material 0

Type 1: PVC/PE Wound Plaster Type 2: Cotton Wound Plaster

Type 3: Elasticity Flexible Cloth Wound Plaster

Type 4: Non-woven Wound Plaster

Type 5: PU Wound Plaster

Type 6: EVA Wound Plaster



White C Skin Blue Green

Transparent

Size 0

⊄ 22mm 72mm x 25mm ⊄25mm 76mm x 25mm 40mm x 10mm 76mm x 38mm 51mm x 45mm(butterfly shape) 38mm x 19mm 60mm x 19mm 76mm x 38mm(H shape) 72mm x 19mm 76mm x 45mm(butterfly shape)




Characteristics:

Highly absorbent and low adherent pad.

· Very hypoallergenic.

Applications: For covering skin scratches and minor wounds.





Type 1: Spun-laced non-woven, with absorbent pad.

Type 2: Cotton cloth, with absorbent pad.

Type 3: Elastic flexible cloth, with absorbent pad.

Type 4: PE/PVC, with absorbent pad



Color 0

White C Skin Pattern



Size 0

4cm x 10cm	8cm x 100cm
6cm x 10cm	4cm x 500cm
8cm x 10cm	6cm x 500cm
4cm x 100cm	8cm x 500cm
6cm x 100cm	



Characteristics:

Highly absorbent and low adherent pad.

Very hypoallergenic.

Applications: For covering skin scratches and minor and larger


WA.

WOUND DRESSING



MED3 Wound Dressing Plaster •







Type 1: Spun-laced Non-woven, With absorbent Pad Type 2: Transparent PU, With absorbent Pad

Type 3: Transparent PU, Without absorbent Pad



O White O Transparent



5cm x 7cm 10cm x 20cm 6cm x 8cm 10cm x 25cm 10cm x 10cm 10cm x 30cm 10cm x 12cm 10cm x 35cm 10cm x 15cm



Characteristics:

- Soft breathable cover.
- Highly absorbent and conformable pad.
- Low allergy adhesive.

Application

Applications:
Ideal for post-operative use, also used in accident & emergency for cuts, laceration and sutured wounds.



MCD41.V. Dressing Plaster O







Type 1: Spun-laced non-woven. Type 2: Transparent PU.



C White C Transparent



6cm x 7cm 6cm x 8cm



Characteristics:

 Soft and comfortable, low allergy adhesive.
 High Permeability, and avoid skin by soaking Applications:


Applications: It can be used for fixing and caring for vein I.V. and all kinds of cannula and infusion lines.

SPENE

WOUND DRESSING


















surgery.


tears. It is also ideal for instances when aesthetic results are of great importance, such as in plastic



WOUND DRESSING











Type 1: Spun-laced non-woven. Type 2: Paper frame transparent film



C White C Transparent




6cm x 9cm (for fixing finger use)


Size & O

Application

- 27cm x 10cm (for finger tip use)
- 7cm x 8.5cm (for fixing arterial and T-long catheter)
- ■10cm x 11.5cm (for fixing multi-lumen catheter,etc) 10cm x 10cm (for drainage use)
- \$12cm x 12cm (for arthrosis wound)
- 611cm x 20cm (for anus & intestine operation)
- Note: We can do other sizes according to customers' requirements.









Hydrocolloid With Film



C Light Yellow



- 44mm x 69mm
- 3 20mm x 60mm
- 28mm x 46mm
- 228mm x 69mm 37mm x 55mm



Size 0

Characteristics:

- Good absorbency for exudation from the wound.
- O Keep the wound moist and accelerate the wound healing, reducing the pain.
- Waterproof, breathing and prevent the wound from bacteria outside.

Applications: Ideal for low or moderate exuded wounds.



MDD1 Absorbent Gauze Roll o



Material O	100% cotton (With or	without X-ray)
<u></u>	8threads/12x8mesh	13threads/19x15mesh
	12threads/19x11mesh	17threads/26x18mesh
Mesh 0	13threads/20x12mesh	20threads/30x20mesh



_			
Town St.		90cm x 100y	90cm x 100m
1	ize O	90cm x 50y	90cm x 50m



Description 0

Characteristics:

Good Whiteness and absorbency.

 Soft, feel comfortable, have good flexibility and micro-pore function.

 Easy to use, and its length is up to you, and also easy to carry.

Applications:

Widely used in the first aid, self-protection and aid in various trainings, sports, fieldwork and family health etc.

MDD2 Non-sterile Gauze Swabs o




Material 0

100% cotton (With or without X-ray, with folded edge or unfolded edge.)

The same of the sa		
6	8threads/12x8mesh	13threads/19x15mesh
	12threads/19x11mesh	17threads/26x18mesh
Mesh 🔾	13threads/20x12mesh	20threads/30x20mesh





10cm x 10cm - 8ply/12ply/16ply



Characteristics:

Good Whiteness and absorbency.

Soft, feel comfortable, have good flexibility and micro-pore function.

Description 0

Applications:

Ideal for wound dressings, wound packing and general wound care.

MDD3 Sterile Gauze Swabs



Material 0

100% cotton (With or without X-ray, with folded edge or unfolded edge.)

	Color 0	C Bleached	
10		12threads/19x11mesh	17threads/26x18mesh
-	Mesh	13threads/20x12mesh	20threads/30x20mesh
		13threads/19x15mesh	

Size 0

5cm x 5cm - 8ply/12ply/16ply 7.5cm x 7.5cm - 8ply/12ply/16ply

10cm x 10cm - 8ply/12ply/16ply



 Good Whiteness and absorbency. O Soft, feel comfortable, tender, have good flexibility and micro-pore function.

Description 0 Applications:

Ideal for wound dressings, wound packing and general wound care. It can be used on wound directly.

SPEND

















general wound care.

Description O O

 Wide range of surgical dressing available to choose from.

Ideal for wound dressings, wound packing and


O Disposable and hygienic.

Applications:

Ideal for wound dressings, wound packing and general wound care. It can be used on wound directly









or blister pack. Applications:

Ideal for trachea fixing and wound care.










petrolatum emulsion to prevent dressing adherence while allowing exudate to easily pass

through the absorbent layer












Characteristics:

 Very smooth and protects the wound, nonadherent. Allows the wound to drain freely into an absorbent secondary dressing.

Applications:

It is used as a primary wound contact layer and the paraffin is present to reduce the adherence of the product to the surface of a granulating wound. Especially for the treatment of burns, ulcers, skin grafts and a variety of traumatic injuries.





Alcohol Swab: Pad saturated with 70% isopropyl alcohol

Alcohol Swab: Pad saturated with 7 with 7 a 10%
 Providone-iodine Pad: Pad saturated with a 10%
 povidone-iodine solution equivalent to 1% available iodine
 Cleaning Wine: One paper with isopropyl alcohol



Alcohol Swab:

30mm x 60mm(2ply), 60mm x 60mm(4ply)

Providone-iodine Pad:

30mm x 60mm(2ply)

Wet Cleaning Wipe:

130mm x 200mm(16ply) 120mm x 150mm(24ply)



Characteristics:

 First aid series, for professional and hospital use.

For external use only.

O Sterile

Applications:

It can be used for cleaning and sterilizing for wounds places and mosquito bites, tableware, children toys, office supplies, etc.





100% cotton (With or without X-ray, with or without rubber ring)



Bleached



20threads/30x20mesh



15cm x 16cm 30cm x 30cm 20cm x 20cm 30cm x 40cm 25cm x 26cm 40cm x 50cm



Characteristics:

Highly soft, good absorbency.

Description 🔾 Applications:

It is mainly used for absorbing blood and exudates during operation and cleaning the wounds.

COTTON



MED1 Absorbent Cotton Wool





Sec.	
1	
	Material 🔾
	Marcella

100% Absorbent Cotton



Bleached



	50g	400g	
	100g	454g	
0	150g	500g	
_	200g	1000g	
	250g		



Characteristics:

 Highly soft, good absorbency.
 Applications:
It is mainly used to absorbing blood and other liquid.



MED2 Zig-zag Cotton Wool o





-	
17 600	
E.	
N. W.	
	Material 0
	Marsia

100% Absorbent Cotton



Bleached Color 0



250g 25g 50g 500g 100g 1000g Size 0 200g



Description 🔾

Characteristics: Soft and hygienically.


Applications: It is suit for cosmetic and medical use.



COTTON








weight, never adhere to mucous membranes.

Applications: It is designed for dentistry surgery.









COTTON







scarred skin appear more natural in color and texture, as well as help to flatten raised scars. It is effective on any part of the body, including the face, and can be worn overnight, while sleeping.















Material 0	10 – 20 g/m2 PP Sp	unbond Non-woven.
Color O	C White	© Green
	18"	21"
Size O	19"	22"
	20"	24"










NON-WOVEN














26

27
















NON-WOVEN



10





Material 0	Type 1: 30/40 g/m² PP S Type 2: 40 - 50 g/m² SM Type 3: 20/25/30 g/m² F	S.
Color	C White C Blue C G	reen C Yellow @ Pink
Color O	White Blue Gi	reen C Yellow C Pink





16









SPENO

29























